

# BMEM Dissociation

Ships at room temperature

Catalog number: 1002

500 mL



## Product Description

BMEM Dissociation is a ready-to-use, highly optimized, chemically defined, animal component-free, enzyme-free dissociation solution for passaging human induced pluripotent stem cells (hiPSCs) in adherent culture.

BMEM Dissociation contains an optimized concentration of dissociation agents and a BMEM-based buffered solution for optimal cell survival during passaging. BMEM Dissociation has been extensively tested for compatibility with weekend-free and no-change feeding schedules.

BMEM Dissociation is recommend for use with hiPSC cultured in BMEM Human Primed Pluripotent (1001) and BMEM Max Pluripotent (1004).

BMEM Dissociation is compatible with a variety of culture matrices, including: Corning® Matrigel® GFR, Gibco™ Geltrex™ GFR, R&D Systems™ Cultrex™ UltiMatrix GFR, SigmaAldrich® ECMatrix™-511 E8 Laminin, Gibco™ Recombinant Vitronectin, Corning® Synthemax II-SC, Acro Biosystems™ Human Laminin 521 Protein

Each lot of BMEM Dissociation is performance tested on hiPSC and tested for sterility.

Component Name	Size	Storage
BMEM Dissociation	500 mL	Store at room temperature

## Instructions for Use

BMEM Dissociation should be stored at room temperature and does not need to be warmed prior to use.

## Directions for Use

- **Step 1**  
When cells reach 70-80% confluency (typically 4 days), aspirate used medium. Add 1 mL of BMEM Dissociation to each well of a 6-well plate (or equivalent) and incubate at room temperature for 1 minute.
- **Step 2**  
After 1 minute aspirate the BMEM Dissociation and leave the plate in the cell culture hood an additional 5 minutes.

**NOTE:** 6 min total is well-optimized and highly recommended.

- **Step 3**  
Using 1 mL of BMEM Human Primed Pluripotent with ROCK inhibitor or BMEM Max Pluripotent with ROCK inhibitor in a P1000 pipette tip, blast the medium against the culture surface to release the cells. Pipette the BMEM/cells up into the P1000 pipette tip (being careful not to make bubbles) and repeat blasting ~5 times until all cells are released.


**CRITICAL:** Avoid pipetting cells too many times as this will decrease cell survival.

- **Step 4**  
Split the cells to the required ratio and transfer to new coated plates

## ROCK Inhibitor

Addition of a ROCK inhibitor, such as 10  $\mu$ M Y27632 or 20 nM Chroman 1 for 24 hours after passage is recommended to provide more reproducible passaging, reduce selective pressure, and allow higher split ratios.

## Representative Images



**Figure 1.** Representative images of passaging BMEM Human Primed Pluripotent-cultured hiPSC using BMEM Dissociation at different time points. **(a)** 80% confluent well before addition of BMEM Dissociation; **(b)** well after the combined 6 minutes of BMEM dissociation and rest; **(c)** immediately after blasting all the cells with BMEM Human Primed Pluripotent with ROCK inhibitor; **(d)** immediately after transferring to the new plate at a 1:25 split ratio; **(e)** 15 min post seeding in the new plate, some cells are already attaching to the bottom of the well. **(f)** 1 h post seeding, the majority of cells attached to the bottom of the well.

## Product Use

For Research Use Only.

Not for diagnostic procedures.